

Q-blot: a highly sensitive immunoblot for detection of *Coxiella burnetii* infections

Anja Scholzen, Katinka D. Mulder, Maarten van der Zanden, Anja Garritsen

Innatoss Laboratories B.V., Oss, The Netherlands



Background

Q fever

- Neglected zoonotic disease with **world-wide prevalence** (1); primarily spread by small ruminants (goat, sheep) and cattle in parturition season
- Caused by the obligate intracellular bacterium *Coxiella burnetii*
- Highly environmentally stable and contagious; airborne transmission** (potential biothreat agent)
- Long-term complications** common (chronic fatigue, persistent infection = chronic Q fever)
- Exposure during **at-risk occupational and domestic activities** (farming, slaughterhouse work, military deployment, rural tourism) (1, 2)



Epidemiology

- Europe:** 824 confirmed notified cases in 2015 (1.9 cases per million); case number in Germany and France increasing since 2013 (3)
- Netherlands:** 2007-2010 outbreak with >4000 notified cases and estimated >40,000 infections (4)
- US:** 2000-2012 0.38 notified cases per million inhabitants per year (5)
- Africa:** studies indicate substantial and underestimated contribution to "fever & infectious endocarditis of unknown origin" (1,6)

Diagnosis/surveillance

- Serology** most common method to test for exposure/infection if suspected (passive surveillance)
- In-house developed tests commonly used (no uniform procedure)
- Indirect immunofluorescence assay (IFA)**
 - Reference method
 - Regarded most sensitive (1)
 - Time-consuming, non-automated, subject to inter-observer variability
- Other commercially available serological tests: ELISA and complement fixation test (CFT)

Objective, Rationale & Approach

Objective

Develop an immunoblot that is more sensitive than IFA and allows fast, automatable, objective detection of exposure to *Coxiella burnetii*

Need for a more sensitive serological Q fever assay

- 2014 field study identified 33,8% of inhabitants in Dutch outbreak village Herpen as IFA-positive (7)
- Parallel assessment by cellular interferon gamma release assay (IGRA, Q-detect™) identified 92% of IFA-positive samples, and 80 % more exposures than IFA (8, 9, unpublished data)
- Question: Can a more sensitive serological assay identify *Coxiella*-specific IgG responses in IGRA-positive, IFA-negative individuals?

Approach & Initial results

- Immunoblot developed by Innatoss Laboratories B.V. using the same proprietary *Coxiella* antigen as Q-detect™ (strain Cb2009-02629 cultured by Wageningen Bioveterinary Research under cell-free conditions)
- First generation Q-blot identified 28/41 (68%) of Q-detect™ IGRA-positive, IFA-negative serum samples from convalescent subclinically/clinically *Coxiella burnetii* infected individuals as sero-positive
- Q-blot was further optimized and tested to address the following question:

Is Q-blot more sensitive than the FocusDX IFA in detecting anti-*Coxiella* IgG?

Q-blot has at least 5-fold higher sensitivity than IFA

Results

- Q-blot signal of hallmark band correlates with phase 2 IgG IFA titers (similar signal across samples at the same dilution as IFA)
- Q-blot shows (lower) **second hallmark band** for all samples with a **phase 1 IgG titer** $\geq 1:128$
- Q-blot is more sensitive than IFA** (still detects anti-*Coxiella* IgG at a 5-fold lower serum dilution than IFA)
- Q-blot is easy to handle;** assay time (manual): ~2 hours 15 minutes; Hands-on time: ~20 minutes




Figure 1: Sensitivity of Q-blot for detection of *Coxiella burnetii* specific IgG with at least 5-fold higher sensitivity than IFA. Heat-killed *Coxiella* antigen was separated using gel electrophoresis followed by Western blotting onto PVDF membranes. Strips were cut using an Accutran Strip Cutter. Incubation with seven different patient sera and detection with polyclonal anti-human-IgG labeled with alkaline phosphatase was performed manually using standard immunoblotting procedures. For each patient sample, 3 strips from the same experiment are shown, with serum samples being diluted (1) (at a fixed dilution factor of 1:51, (2) the dilution factor of the phase 2 IgG IFA titer or (3) 5-fold higher than the phase 2 IgG IFA titer of the respective sample. Phase 1 IgG IFA titers were determined at the Jeroen Bosch Hospital, Den Bosch, NL using the FocusDX IFA. IFA = Indirect immunofluorescence assay; QFS = Q fever fatigue syndrome

Ongoing work

- Conversion** from research-only-use in house assay to **commercial product**
 - Development of **technically stable product**
 - Test **lot-to-lot variation** using different *Coxiella* antigen batches
 - Establishment of **automated** protocol
- Clinical validation against IFA and Q-detect™ IGRA** to determine sensitivity and specificity using an in-house collection of sera from
 - Convalescent individuals with past acute or asymptomatic Q fever
 - Chronic Q fever patients
 - Q fever fatigue syndrome patients

Outlook

Clinical utility in Q fever surveillance

- Sensitive, fast and easy to use assay
- Suitable for remote locations where implementation of the Q-detect™ IGRA is less feasible
- Suitable for (semi) automated screening of large serum collections

Looking for collaborations and customers!



References

- (1) Eldin et al. From Q Fever to *Coxiella burnetii* Infection: a Paradigm Change. Clinical Microbiology Reviews. Jan 2017.
- (2) Farris et al. Seroconversions for *Coxiella* and Rickettsial Pathogens among US Marines Deployed to Afghanistan, 2001-2010. Emerging Infectious Diseases. Aug 2016.
- (3) European Centre for Disease Prevention and Control. ECDC annual epidemiological report for 2015 – Q-fever. Stockholm Dec 2017.
- (4) Kampichore et al. Screening for *Coxiella burnetii* seroprevalence in chronic Q fever high-risk groups reveals the magnitude of the Dutch Q fever outbreak. Epidemiology and Infection. April 2013.
- (5) Dahlgren et al. Q Fever in the United States: Summary of Case Reports from Two National Surveillance Systems, 2000–2012. American Journal of Tropical Medicine and Hygiene. Feb 2015.
- (6) Vanderburg et al. Epidemiology of *Coxiella burnetii* Infection in Africa: A OneHealth Systematic Review. PLoS Neglected Tropical Diseases. April 2014.
- (7) Morroy et al. Population Screening for Chronic Q Fever Seven Years after a Major Outbreak. PLoS One. July 2015.
- (8) Schoffelen et al. Specific interferon gamma detection for the diagnosis of previous Q fever. Clinical Infectious Diseases. Jun 2013.
- (9) Garritsen et al. Experiences with detection of *Coxiella* exposure using Q-detect™, an IGRA for Q fever. ESCMID conference. 2017.